

Diagnopal Dyanabioassay Products LTD

HRP- Conjugates Stabilizer (BSA Free)

DPHRPS-02 (100ml), DPHRPS-02 (1000ml),

Introduction:

Diagnopal HRP-conjugate stabilizer (**DPHRPS-02**) is proper choice to make ready to use HRP-Ab, HRP-Protein, and HRP-peptide conjugates for immunoassay. It is recommended to use for ELISA and chemiluminescence. Standard recovery and linearity for immunoassay tests, blocking of proteins and antibodies Non-Specific Bindings and long-term stability of conjugates is provided by DPHRPS-02 for your tests. It makes more sensitive and specific Immunoassay results. Due to its unique production method, DPHRPS-02 will help you to produce more reproducible products and more repeatable test results in your commercial Kits.

Specification of “DPHRPS-02 stabilizer”


Stability and storage	The product should be kept at 2-8 °C. In refrigerated condition it is stable for 3 years from date of manufacture.
Specification	Protein: BSA Free
	Product Buffer: NONE
	pH: 6.3-7.3
	Preservative: Gentamycin, Methyl Isothiazolone
	Acceptable Diluents: Do not dilute
	Dry weight(gr/mL): 0.052-0.076
Shipment Temperature: RT	
Notes	Keep product away from freezing, Direct sun light and heat. 3 months of storage of un-opened products in Dark and RT will not affect the product quality or efficacy more than 0.05%.
Available forms	100 ML and 1Litter standard vials and requested volume.

Diagnopal Dyanabioassay Products LTD

Recommendation for use:

DPHRPS-02 is designed to be used as stabilizer for HRP conjugated proteins, Abs and Ags.

Stability of Conjugates in DPHRPS-02



Related Products:

[DPHRPS-01](#)

[DPHRPS-02](#)