

# Diagnopal Dynabioassay Products LTD

## HRP- Conjugates Stabilizer (protein free)

DPHRPS-01 (100ml), DPHRPS-01 (1000ml),

### Introduction:

Diagnopal HRP-conjugate stabilizer (**DPHRPS-01-**) is proper choice to make ready to use HRP-Ab, HRP-protein and HRP-peptide conjugates for immunoassay. It is recommended to use for ELISA, CLIA. Standard recovery and linearity for immunoassay tests. Diagnopal HRP-conjugate stabilizer (**DPHRPS-01**) block proteins interferences and antibodies, Non-Specific Binding and prepares long term stability of conjugates for your tests. It makes more sensitive and specific Immunoassay results and due to its unique production method, DPHRPS-01 will help you to produce more reproducible products and more repeatable test results in your commercial Kits.

### Specification of “DPHRPS-01 stabilizer”


<b>Stability and storage</b>	The product should be kept at 2-8 °C. In refrigerated condition it is stable for 3 years from date of manufacture.
<b>Specification</b>	Protein: NONE
	Product Buffer: NONE
	pH: 6.3-7.3
	Preservative: Gentamycin, Methyl Isothiazolone
	Acceptable Diluents: Do not dilute
	Dry weight(gr/mL): 0.052-0.074
Shipment Temperature: RT	
<b>Notes</b>	Keep product away from freezing, Direct sun light and heat. 3 months of storage of un-opened products in Dark and RT will not affect the product quality or efficacy more than 0.05%.
<b>Available forms</b>	100 ML and 1Litter standard vials and requested volume.

# Diagnopal Dynabioassay Products LTD

## Recommendation for use:

DPHRPS-01 is designed to use as stabilizer for HRP conjugated proteins, Abs and Ags.

## Stability of Conjugates in DPHRPS-01



Related Products:

DPHRPS  
DPHRPS-01