

Diagnopal Dynabioassay Products LTD

Daignopal Hapten/Hapten-HRP Conjugate Stabilizer

General Hapten, hapten/HRP- Conjugated Stabilizer

DPHHRPG (100ml), DPHHRPG 1000 (1000ml)

Introduction:

Diagnopal Hapten/ Hapten-HRP-conjugate stabilizer (**DPHHRPG**) is proper choice to make ready to use Hapten, Hapten/HRP-conjugate s for immunoassay. It is recommended for ELISA, RIA and chemiluminescence and or making ready to use controls and standards for immunoassay test for long-term storing in refrigerator. Standard recovery and linearity for immunoassay tests, Non-Specific Binding and long-term stability of conjugates is provided by DPHHRPG for your tests and make more sensitive and specific Immunoassay results. Due to its unique production method, DPHHRPG will help you to produce more reproducible products and more repeatable test results in your commercial Kits or research. This product has been tested for T3, T4, FT3, FT4, T-Uptake, Progesterone, Testosterone, DHEA, fluoroquinolone antibiotics, Aminoglycosides, HCG, FSH, TSH and some other Haptens.

Specification of “DPHHRPG stabilizer”


Stability and storage	The product should be kept at 2-8 °C. In refrigerated condition it is stable for 3 years from date of manufacture.
Specification	Protein: BSA
	Product Buffer: None
	pH: 6.4-7.3
	Preservative: Gentamycin, Methyl Isothiazolone
	Acceptable Diluents: Do not dilute
	Dry weight(gr/mL): 0.050-0.075
	Shipment Temperature: RT
Notes	Keep product away from freezing, Direct sun light and heat. 4 months of storage of un-opened products in Dark and RT will not affect the product quality or efficacy more than 0.05%.

Diagnopal Dynabioassay Products LTD

Recommendation for use:

DPHHRPG is designed to use as stabilizer for Hapten, Hapten-HRP conjugates such as T3-HRP, T4-HRP, Drugs/Drug-HRP and etc. For more products please see our web page.

Stability of Conjugates in DPHHRPG



Related Products: Please see our web page.

This product is for research and/or further manufacture use only and is not for human therapeutic or consume purpose. Sales without any seller's warranty or representation. No claims for replacement or refund of purchase shall be allowed. All claims must be made within 30 following date of delivery.